

Chart1: Laboratory Approach to Etiological diagnosis of ADD

Fresh stool or Rectal swab/ Fecal swab in Cary Blair

Macroscopic Examination
(of Stool)

Examine and record consistency
(example: mucous, blood, liquid)

Microscopic examination

Chart 2

Culture

Direct Plating
MAC, XLD or HE, and
TCBS

APW S/c after 4- 8hrs
SF S/c after 16-24 hrs

Enrichment
(APW, SF broth)
At 36+/-1°C

16-18hrs at 36°C +/-1°C

TCBS

Yellow, Shiny colonies
(1-4 mm)

Chart 3

MAC

Colorless (NLF), smooth
colonies (1-4 mm)

Chart 4

XLD/HE

XLD: Red colonies with or
without black centre (1-3)
mm
HE: Colourless colonies
with or without black centre
(1-3)mm

MAC= MacConkey agar, XLD= Xylose Lysine Deoxychlate agar, TCBS= Thiosulphate Citrate Bile salt
Sucrose agar, APW= Alkaline Peptone Water, SF= Selenite F broth

Chart 2: Microscopic Examination of Stool/Rectal Swab/Fecal Swab


Preparation of Saline and Iodine mounts:

- Place 1 drop of 0.85% NaCl on left side of a slide (1- x 3-in) and 1 drop of iodine (5 times diluted Lugol's iodine) on right side of slide
- Take small amount of fecal specimen (using a swab/loop) and thoroughly mix
- Place a 22-mm coverslip on each suspension




Method of examining direct wet film preparation with a 10x objective:

- Scan both the suspensions with 10x objective.
- The entire coverslip area should be examined for motile parasites under low power
- If something suspicious seen-- 40x objective for more detailed study
- At least 1/3 of coverslip should be examined under 40x even if nothing suspicious seen




NaCl: Sodium Chloride; Lugols Iodine: I_2 : KI: H_2O = 5 gm:10 gm:100ml

Chart3: Identification of *Vibrio cholerae*



*Antibiotic discs required: Ampicillin (10µg), Tetracycline (30µg), Co-trimoxazole (1.25/23.75µg), Ciprofloxacin (5µg)

Chart 4: Identification of Salmonella & Shigella



XLD: Pink/Red colonies ± black center (2-4 mm)
HE: Clear colonies ± black center (1-4mm)

Subculture on Nonselective agar (MHA or NA)

Oxidase test

TSI, LIA, MIO, Citrate, Urea

Antimicrobial
Susceptibility Testing*

TSI:K/A, g++
LIA: LDC +
MIO: +/-/+
Citrate:+
Urea: -

TSI: K/A tr
LIA: LDC +
MIO: +/-/-
Citrate:-
Urea: -

TSI: K/Ag
LIA: LDC -
MIO: +/-/+
Citrate:-
Urea: -

TSI: K/A
LIA: LDC -
MIO: -/v/-
Citrate:-
Urea: -

TSI: K/A
LIA: LDC -
MIO: -/-/+
Citrate:-
Urea: -

Consistent with
Salmonella
species

Consistent with
Salmonella
serov. Typhi

Consistent with
Salmonella ser.
Paratyphi A

Consistent with
Shigella species

Consistent with
Shigella sonnei

*Antibiotic discs required: Ampicillin (10µg), Co-trimoxazole(1.25/23.75µg), Chloramphenicol (30µg), Ceftriaxone (30µg),Pefloxacin (5µg),Ciprofloxacin (5µg)

Chart5 : Disc Diffusion Method for Susceptibility Testing

Isolated colonies



Saline Suspension



Adjust Turbidity to 0.5 Mc Farland



Inoculate Mueller Hinton Agar
Inoculate a quality control
plate in the similar manner
using
E. Coli 25922





Place discs on Plate



Incubate 16-18 hr,
 $36 \pm 1^{\circ}\text{C}$

Measure zone sizes,
Interpret



Antibiotic discs required:
Salmonella & *Shigella*: Ampicillin (10 μg),
Co-trimoxazole(1.25/23.75 μg),
Chloramphenicol (30 μg), Ceftriaxone
(30 μg),Pefloxacin (5 μg),Ciprofloxacin
(5 μg)

Vibrio cholerae: Ampicillin (10 μg),
Tetracycline (30 μg), Co-trimoxazole
(1.25/23.75 μg), Ciprofloxacin (5 μg)

Antibiotic disc placement:
Discs must be placed so that unacceptable
overlapping of zones is avoided.

The maximum number of discs will vary
based on the antimicrobial and organism
(some agents consistently give larger zones
with susceptible isolates).

No more than six (6) disks should be
placed on a 90mm plate and no more than
twelve (12) disks should be placed onto a
150 mm plate

Triple Sugar Iron (TSI) agar

- Principle:

To determine following characteristics in a GNB

- Fermentation of Glucose, Sucrose & Lactose
- Gas production during glucose fermentation
- H₂S production


- Medium:

- Triple Sugar Iron agar containing Glucose, Sucrose, Lactose (1:10:10), Phenol red indicator, Iron salt, Sodium thiosulfate (H₂S source)

- Procedure:

- Inoculate the agar heavily (Stab and Streak)
- Incubate at 37°C, examine after 8-12 hrs, & 18-24 hrs
- If screw capped tubes are used, the caps must be loose.

Triple Sugar Iron (TSI) agar



Reactions: After at 18-24 hours:

- A: Uninoculated tube
- B: K/A Tr(ace) H₂S: glucose fermented, lactose sucrose not fermented. *Salmonella* Typhi
- C: K/Ag++: Glucose fermentation but no sucrose or lactose fermentation. Non-typhoidal *Salmonella* spp.
- D: K/A: Glucose fermented, no lactose or sucrose fermented, no gas produced (*Shigella flexneri*).
- E: A/A: Glucose and lactose and/or sucrose fermented gas produced. *E. coli*.
- F: No fermentation, no gas production. *Pseudomonas aeruginosa*.

Citrate Utilization test

Principle:

- Test the ability to utilize citrate as a sole carbon source

Medium:

- Simmon's citrate medium (solid medium with Bromothymol blue indicator)

Procedure:


- Light inoculum (Streak)
- Incubation at 37°C, up to 48 hrs

Result Interpretation:

- Positive: Growth & blue color
- Negative: No growth, green color (no change)

Controls:

- Positive: *Citrobacter* spp
- Negative: *E. coli*




Urease Test

- **Principle:**
 - Test for the presence of enzyme Urease that breaks urea into NH_3 , thus increasing the pH
- **Medium:**
 - Cristensen's urea agar, Phenol red indicator
- **Procedure:**
 - Inoculate the agar heavily
 - Incubate at $36^\circ \pm 1^\circ \text{C}$, examine after 4 hours, after overnight incubation up to 4 days
- **Result Interpretation:**
 - Pink colour -----Positive
 - No coloration-----Negative
- **Controls:**
 - Positive: *Proteus spp.*
 - Negative: *E. coli*



Motility Indole Ornithine (MIO) agar



A

B

C

D

E

F

Reactions

A: Uninoculated MIO tube

B: Indole negative (*Shigella flexneri*)

C: Indole positive (*E. coli*)

D: Motile/Ornithine Positive (*Salmonella* serovar Newport)

E: Nonmotile/Ornithine Negative (*Shigella flexneri*)

F: Nonmotile/Ornithine Positive (*Shigella sonnei*)

Motility Indole Ornithine (MIO) agar

Principle:

- To demonstrate the motility, ornithine decarboxylase activity and Indole production.

Procedure:

- Inoculate the agar with a single stab into agar and stop approx.1cm from the bottom of the tube.
- Incubate at 37°C, examine after overnight (18-24 hrs) incubation .

Result Interpretation::

Motility

- Positive: Turbid.
- Negative: Clear

Motility Indole Ornithine (MIO) agar Cont.....

Ornithine Decarboxylase

- Positive: middle of tube turns light purple colour
- Negative: Yellow


Indole: (Kovac's reagent)

- Positive: Pink Red
- Negative: No Color Change

Controls:

- Positive: *E.coli* spp.
- Negative: *Shigella flexneri*

Lysine Iron Agar (LIA)



Reactions

A: Uninoculated LIA

B: LDC negative /lysine deaminase negative/ H_2S negative (*Citrobacter freundii*)

C: LDC positive /lysine deaminase negative/ H_2S negative (*Salmonella* ser. Typhi)

D: LDC negative /lysine deaminase positive/ H_2S negative (*Proteus mirabilis*)

E: LDC positive /lysine deaminase negative/ H_2S positive (*Salmonella* ser. Newport)

Lysine Iron Agar (LIA)

Principle:

To demonstrate the Lysine decarboxylation, deamination and H₂S production by enteric organisms.

Procedure:

- Inoculate the agar with a single stab into agar and stop approx.1cm from the bottom of the tube.
- Incubate at 37°C, examine after overnight (18-24 hrs)incubation .

Controls:

- *Proteus mirabilis*: Slant Red; Butt Yellow; H₂S neg
- *Salmonella* Typhimurium: Slant Purple; Butt Purple; H₂S pos
- *Shigella flexneri*: Slant Purple; Butt Yellow; H₂S neg

Lysine Iron Agar (LIA) Cont....

Result Interpretation:

H₂S Production

- Positive: black colour along the streak or throughout medium
- Negative: No Colour change

Lysine Decarboxylation (detected in butt of the tube)

- Positive: Purple
- Negative: Yellow

Lysine Deamination (detected on agar slant)

Positive: Red

Negative: Purple

Salmonella serovar Typhi



A

B

C

D

E

F

A) TSI: Alkaline slant / Acid Butt / Trace H₂S / No Gas (K / A TR)

B) Urea: Negative

C) LIA: Lysine Decarboxylase Positive

D) Citrate: Negative

E) MIO: Motile / Ornithine Negative

F) MIO w/ indol reagent: Indol negative

***Salmonella* serovar Newport (typical of most nontyphoidal *Salmonellae*)**



A

B

C

D

E

F

A) TSI: (K/Ag++) Alkaline slant / Acid Butt / H₂S Positive / Gas

B) Urea: Negative

C) LIA: Lysine Decarboxylase Positive

D) Citrate: Positive

E) MIO: Motile / Ornithine Negative

F) MIO w/ indole reagent: Indole negative

Salmonella serovar Paratyphi A



A

B

C

D

E

F

A) TSI: Alkaline slant / Acid Butt / No H₂S / Gas (K / A g)

B) Urea: Negative

C) LIA: Lysine Decarboxylase Negative

D) Citrate: Negative

E) MIO: Motile / Ornithine Positive

F) MIO w/ indole reagent: Indole negative

Shigella sonnei



A

B

C

D

E

F

A) TSI: Alkaline slant / Acid Butt / No H₂S / No Gas (K / A)

B) Urea: Negative

C) LIA: Lysine Decarboxylase Negative

D) Citrate: Negative

E) MIO: Nonmotile / Ornithine Positive

F) MIO w/ indole reagent: Indole negative

Shigella flexneri



A

B

C

D

E

F

A) TSI: Alkaline slant / Acid Butt / No H₂S / No Gas (K / A)

B) Urea: Negative

C) LIA: Lysine Decarboxylase Negative

D) Citrate: Negative

E) MIO: Nonmotile / Ornithine Negative

F) MIO w/ indol reagent: Indol negative

Escherichia coli



A

B

C

D

E

F

A) TSI: Acid slant / Acid Butt / No H₂S / Copious Gas (A / A g)

B) Urea: Negative

C) LIA: Lysine Decarboxylase Positive

D) Citrate: Negative

E) MIO: Motile / Ornithine Positive

F) MIO w/ indol reagent: Indol positive

References

- Laboratory Protocol: “Isolation of Salmonella and Shigella from Faecal Specimens”
<http://www.antimicrobialresistance.dk/data/images/protocols/gfn%20stool%20culture-nov2010.pdf>
- Laboratory Protocol: “Biochemical Identification of Salmonella and Shigella Using an Abbreviated Panel of Tests”
http://www.antimicrobialresistance.dk/data/images/protocols/gfn_biochem_final.pdf
- Laboratory Protocol: “Susceptibility testing of Enterobacteriaceae using disk diffusion”
http://www.antimicrobialresistance.dk/data/images/protocols/sop_dd_2010gfnlab002_v2.pdf